

# VISIOCAM 4

communicate with sharp images



WITH VISIOCAM 4 CAMERA, COMMUNICATION, MOTIVATION AND LOYALTY DEVELOPMENT WITH YOUR PATIENTS HAVE NEVER BEEN SO EASY. YOU CAN VISUALISE INSTANTANEOUSLY YOUR PICTURES ON YOUR COMPUTER. PICTURES ARE PURE AND COLOURS ARE SURPRISING OF REALITY. YOUR PATIENT WILL UNDERSTAND IMMEDIATELY THE INTEREST OF HIS TREATMENT !



### 6 LED LIGHTING

Thanks to its 6 integrated LED, and its defogging features, the Visiocalm 3 has never been so reliable.

### ULTRA-LIGHT HANDPIECE

You will be impressed by the lightness of its handpiece (38 g).



### USB DIRECT

Connect your Visiocalm 4 directly with a USB cable. Now, you do not need to connect your camera to a power supply.



### SHARP IMAGES FROM SMILE TO MACRO

Thanks to its adjustment button, you can control the sharpness of your image easily. The image is always clear thanks to a large depth of field (5 to 30mm) and high-quality optics (covering a field of 80°). You benefit from a high-quality image for intraoral pictures (macro) or smile pictures.

#### Specification

CCD	1/4" from Sony	LENGTH OF CABLE	3 m / Option: 5 m (USB 2.0)
RESOLUTION	768 x 494 px	DIMENSION	(W) 22 x (D) 22 x (H) 190 mm
LIGHT SOURCE	6 White LEDs	HAND PIECE WEIGHT	36 g
FOCUS RANGE	Smile to Close up	RATING	5 V, 500 mA
ANGLE OF VIEW	80°	POWER SUPPLY	No separate power supply required (via USB2.0 port)



## VIVEZ VOS IMAGES SUR VISIODENT VISIODENT

Visualise with the Imaging software, file your pictures easily on the patient's dossier and give a new dimension to your communication. With retouch tools, it is easier to explain treatment options to your patient. Video-diagnostic has never been so simple!



[www.visiodent.com](http://www.visiodent.com) EMAIL / [contact@visiodent.com](mailto:contact@visiodent.com)

**VISIODENT** IMAGINE & CREATE